

Restriction Map of pTet-Off Vector. Unique restriction sites are in bold.

Description

pTet-Off expresses the tet-responsive transcriptional activator (tTA) from the strong immediate early promoter of cytomegalovirus (P_{CMV}). tTA is a fusion of amino acids 1–207 of the tet repressor (TetR) and the negatively charged C-terminal activation domain (130 amino acids) of the VP16 protein of herpes simplex virus. pTet-Off was originally described as pUHD15-1neo by Resnitzky *et al.* (1994) and was created by insertion of a neomycin resistance gene into pUHD15-1 (Gossen & Bujard, 1992). pTet-Off can be distinguished from pTet-On by digestion with *Hind* III.

Use

The pTet-Off Vector is used to develop stable Tet-Off cell lines. After a vector that contains a gene under the control of a tet-responsive element (TRE) is transfected into a Tet-Off cell line, the tTA binds to the TRE, thus activating transcription in the absence of tetracycline (Tc), or its derivative doxycycline (Dox). As Tc or Dox is added to the culture medium, transcription from the TRE is turned off in a highly dose-dependent manner. More information on TRE-containing vectors and protocols describing the construction of Tet-Off cell lines can be found in the Tet Systems User Manual (PT3001-1).



United States/Canada 800.662.2566 Asia Pacific +1.650.919.7300 Europe +33.(0)1.3904.6880 Japan +81.(0)77.543.6116

Clontech Laboratories, Inc. ATakara Bio Company 1290 Terra Bella Ave. Mountain View, CA 94043 Technical Support (US) E-mail: tech@clontech.com www.clontech.com

(PR23072; published 19 March 2002)

Location of Features

- Fragment containing *P*_{CMV} : 86–673
- Tetracycline-responsive transcriptional activator (tTA): 774–1781
- Col E1 origin of replication: 2604-3247
- Ampicillin resistance gene: β-lactamase coding sequences: 4255–3395
- Fragment containing the SV40 poly A signal: 1797-2254
- Neomycin/kanamycin resistance gene: 6462-5668
- SV40 promoter (P_{SV40}) controlling expression of neomycin/kanamycin resistance gene: 7125–6782.

Propagation in *E. coli*

- Suitable host strains: DH5 α and other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (50 µg/ml) on E. coli hosts.
- E. coli replication origin: Col E1

References

- 1. Tet Expression Systems and Cell Lines (July 1996) Clontechniques XI(3):2-5.
- 2. Gossen, M. & Bujard, H. (1992) Proc. Natl. Acad. Sci. USA 89:5547-5551.
- 3. Gossen, M., et al. (1995) Science **268**:1766–1769.
- 4. Resnitzky, D., et al. (1994) Mol. Cell. Biol. 14:1669–1679.

Note:

The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced.

Notice to Purchaser

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes, nor is it intended for human use. Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of Clontech Laboratories, Inc.

Use of the Tetracycline controllable expression systems (the "Tet Technology") is covered by a series of patents including U.S. patents No. 5,464,758 and No. 5,814,618, which are proprietary to TET Systems Holding GmbH & Co. KG. Academic research institutions are granted an automatic license with the purchase of this product to use the Tet Technology only for internal, academic research purposes, which license specifically excludes the right to sell, or otherwise transfer, the Tet Technology or its component parts to third parties. Notwithstanding the above, academic and not-for profit research institutions who's research using the Tet Technology is sponsored by for profit organizations, which shall receive ownership to all data and results stemming from the sponsored research, shall need a commercial license agreement from IP Merchandisers in order to use the Tet Technology. In accepting this license, all users acknowledge that the Tet Technology is experimental in nature. TET Systems Holding GmbH & Co. KG makes no warranties, express or implied or of any kind, and hereby disclaims any warranties, representations, or guarantees of any kind as to the Tet Technology, patents, or products. All others are invited to request a license from TET Systems Holding GmbH & Co. KG prior to purchasing these reagents or using them for any purpose. Clontech is required by its licensing agreement to submit a report of all purchasers of the Tet-controllable expression system to IP Merchandisers, Inc. For license information, please contact:

Hans Peter Kneubuehl TET Systems Holding GmbH & Co. KG Im Neuenheimer Feld 582 69120 Heidelberg Germany Tel +49 6221 588 04 00 Fax +49 6221 588 04 04 eMail: kneubuehl@tet-systems.de or use our electronic licensing request form via http://www.tetsystems.com/main_inquiry.htm

Clontech, Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc. Clontech is a Takara Bio Company. ©2005