

Use

RNAi-Ready pSIREN-RetroQ is used for targeted gene silencing when an oligonucleotide encoding an appropriate shRNA is inserted. To construct recombinant pSIREN-RetroQ, first design, generate, and anneal complementary shRNA oligonucleotides using the protocols in the Knockout RNAi Systems User Manual (PT3739-1). The annealed oligonucleotide should contain 5'-*Bam*H I and 3'-*Eco*R I sites. Then ligate the annealed oligonucleotide into RNAi-Ready pSIREN-RetroQ.

Your pSIREN-RetroQ construct can be transfected as a plasmid expression vector to screen for functional shRNA oligonucleotides. For gene silencing experiments using viral delivery, transfect the pSIREN-RetroQ construct into a packaging cell line (such as the RetroPack™ PT67 Cell Line (Cat. No. 631510), AmphiPack™-293 Cell Line, EcoPack2™-293 Cell Line, or Pantropic System). RNA from the vector is packaged into infectious retroviral particles. These infectious particles are replication-incompetent since pSIREN-RetroQ lacks structural genes (*gag*, *pol*, and *env*) necessary for particle formation and replication; however, these genes are stably integrated as part of the packaging cell genome. These retroviral particles can infect a wide range of target cells and transmit the shRNA but cannot replicate within these cells due to the absence of viral structural genes. The separate introduction and integration of the structural genes into the packaging cell line minimizes the chances of producing replication-competent virus due to recombination events during cell proliferation.

Location of Features

- PGK promoter (P_{PGK}): 2–510
- Puromycin resistance (Puro^r): 531–1130
- 3' MoMuLV LTR(deletion in U3): 1311–1736
Poly A region: 1562–1577
- SV40 promoter: 2016–2302
- SV40 ori: 2237–2302
- Col E1 ori (Site of replication initiation): 3221–2622
- Ampicillin resistance gene (β -lactamase): 4243–3383
Start codon (ATG): 4243–4241 stop codon (TAA): 3385–3383
- 5' LTR (CMV/MSV): 4606–5333
Cytomegalovirus (CMV)/ mouse sarcoma virus (MSV) hybrid promoter: 4606–5116
R region: 5189–5259
U5 region: 5261–5333
- Ψ^+ (extended packaging signal): 5363–6172
- Human U6 promoter (P_{U6}): 6188–6445

Sequencing Primer Location

- U6 Forward Sequencing Primer: 6188–6206
5'-GGGCAGGAAGAGGGCCTAT-3'

Propagation in *E. coli*

- Suitable host strains: DH5 α , DH10B, and other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (100 μ g/ml) in *E. coli* hosts.
- *E. coli* replication origin: Col E1
- Copy number: low

References

1. Kinsella, T. M. & Nolan G. P. (1996) *Hum. Gene Ther.* **7**:1405–1413.
2. Ory, D. S., Neugeboren, B. A. & Mulligan, R. C. (1996) *Proc. Nat. Acad. Sci. USA* **93**:11400–11406.
3. Pear, W. S., Nolan, G. P., Scott, M. L. & Baltimore, D. (1993) *Proc. Natl. Acad. Sci. USA* **90**(18):8392–8396.
4. Yang, S., Delgado, R., King, S. R., Woffendin, C., Barker, C. S., Yang, Z. Y., Xu, L., Nolan, G. P. & Nabel, G. J. (1999) *Hum. Gene Ther.* **10**:123–132.
5. Emerman, M. & Temin, H. M. (1984) *Cell* **39**:449–467.

Note: The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced.

The viral supernatants produced by this retroviral vector could, depending on your cloned insert, contain potentially hazardous recombinant virus. Due caution must be exercised in the production and handling of recombinant retrovirus. Appropriate NIH, regional, and institutional guidelines apply.

Notice to Purchaser

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of Clontech Laboratories, Inc.

Clontech, Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc. Clontech is a Takara Bio Company. ©2006